



IMPORTANT!

How to Make Sure You Get the Right Sound Level for Your Needs

1. Start by setting the **Tone** control in the middle between the “Low” and “High” settings.



2. Make sure the **Volume** control is in the middle between the “Low” and “High” settings.



3. Place a call.

4. Press the **VOL** button to turn on amplification.



5. As you listen to the person on the other end of the line, experiment by **slowly** increasing the **Volume** setting. Also, try moving the **Tone** control slightly toward the “Low” or “High” settings.

You need to set this only once. Every time you press the **VOL** button, the CrystalTone will use these settings.

NOTE: Some hearing aids may be sensitive when high-frequency sounds are amplified. If you experience a “crackling” sound, please try the following:

- Move the **Tone** control slightly toward “Low” until the crackling stops.
- Move the **Volume** control slightly toward “Low” until the crackling stops.
- Change the **Crystal Tone** switch to the “B” setting to reduce high-frequency amplification. Then, repeat the steps above. (This switch is located on the bottom of your phone.)